# DAYVIA®



LUMINOTHERAPIE

---

## Lampe de luminothérapie DAYVIA **Sundesk / 02**

- NOTICE D'UTILISATION P.2
- BEDIENUNGSANLEITUNG P.9
- GEBRUIKSAANWIJZING P.15
- USER'S GUIDE P.21

## INTRODUCTION

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
# USER'S GUIDE

# INTRODUCTION

Congratulations for purchasing this DAYVIA light therapy lamp.

Used in correct conditions, it will enable you to regain the light that you are lacking.

For certain people, the lack of light can cause various symptoms:

- Irritability and sullen mood.
- Significant increase in appetite.
- Decrease of activity and the efficiency in work due to a lack of energy.
- Sleeping disturbances, chronic tiredness.
- Sexual disturbances (loss of libido).

Many international studies show that light therapy contributes to regain energy, form, ability to concentrate and good mood, the same that you have whilst sunbathing.

The light therapy lamp also contributes to calming the sleeping rhythm for people who do night shifts, people who suffer from delayed sleep phase syndrome (falling asleep late and difficulties waking up in the morning) or advanced sleep phase syndrome (waking up too early) or people who travel and go cross many time zones (jet-lag).

The effect of the light will depend on its intensity but especially on its spectral curve (color). Indeed, certain receptors located in our eyes are very sensitive to a light comprising certain wavelengths just below and above 500 nanometers (blue in particular and / or green), they have a major role on our biological clock. Thus, a blue monochromatic light can be as effective as a much more intense white light.

Your device therefore proposes 2 possibilities as effective as each other:

- Intense white light
- A monochromatic blue light

Please read this user's guide attentively in order to benefit from all the advantages of your DAYVIA lamp.

# SAFETY MEASURES

-The device must be installed and operated in conformity with the information provided on the accompanying documents.

-Do not use the DAYVIA lamp in a damp environment (for example in a bathroom, next to a bathtub or shower cubicle, outside...).

-Keep the power cord away from heated surfaces.

-Only use the supplied DAYVIA power supply unit with your light therapy lamp. This power supply unit is specifically made for your lamp, it must not be used with any other equipment.

-The plug serves as the power selector. The device must be used with an accessible plug.

-This lamp must not be used by children, or physically or mentally disabled people.

-Any modifications to this lamp are forbidden.

-The device must not be used if the cord and/or the power supply unit are apparently damaged or defective.

-The lamp requires special precautions, it can be installed in all rooms including domestic rooms and away from HF sources where the intensity of electromagnetic disturbances is high.

-Avoid using this unit next to other equipment or stacked with other equipment because this may cause malfunction. If this is necessary, this appliance and other appliances should be observed to check for normal operation.  
-Do not use RF portable communication devices (including peripherals such as antenna cables and external antennas) closer than 30 cm (12 inches) to any part of the DAYVIA lamp. Otherwise, the performance of these devices may be impaired.

-You are undergoing a treatment with psychotropic or photosensitising medication

If you encounter any problems (for example headaches or eye tiredness) after the use of the light therapy lamp, please tell your general practitioner. Often, all you have to do is start the treatment with shorter sessions (10 minutes) and gradually increase them each day until the recommended duration.

## PRECAUTIONS

**The light therapy lamp must be used in an environment lit by another source of light,** the lamp is complementary to normal light (electric or natural light).

Contrary to natural daylight or to other lamps on the market, this light therapy medical device diffuses no UV rays.

Studies have shown that light therapy can be used in safety (for example your light therapy lamp diffuses a monochromatic blue or white intense light up to 10000 lux for the white, whereas on a beach in the summer, the light can reach up to 100000lux with on the top harmful rays such as UVs of which your lamp has none).

Nevertheless, before starting your treatment, please consult your general practitioner if :

- You suffer or have suffered from depression
- You are hypersensitive to light
- You suffer from an eye illness
- You suffer from a mental disorder/condition

## TURNING ON/OFF

The Dayvia Sundesk lamp is carefully controlled before packaging and expedition.

You will find the user's guide, the power supply unit and your lamp in the package.

We recommend you keep the packaging of a return is necessary, then it must be discarded in conformity with your country's regulations.

Here are the different steps to turn on your DAYVIA light therapy lamp:

- 1) Insert the male connector coming out of the lamp into the female connector situated on the supplied power supply unit.
- 2) Check that the tension indicated on the power supply unit corresponds to the local sector tension, then plug in the device.
- 3) Turn on the lamp, in a lit environment, using the on/off switch located on the top right-hand side of the lamp.

You can then position the illuminating part to your face during the light therapy session and then modulate the intensity with the - and + keys located next to the switch. The first push on the switch triggers an intense white light which is

very effective from a therapeutic point of view. When you reduce the intensity with the - button, the light will become warmer for office lighting after your light therapy session. A second push on the on / off triggers a low intensity monochromatic blue light as effective as the intense white light.

4) When finished using the device, turn it off by simply pressing the on/off switch.

5) Unplug the device from the socket.

## USER'S RECOMMENDATIONS

Place the lamp between 30 and 70cm from your eyes in the field of vision positioning the illuminating part towards the face during the session, it's not necessary to continually look at the lamp, pursue your activities normally (breakfast, reading, office work, computer, make-up...)

For a pleasant use, without a dazzling effect, make sure that the room is lit (especially for the first days of use) to minimize the contrast between the lighting of the lamp and that of the environment.


It is advised to start the light therapy sessions as soon as you feel a lack of light, before the first symptoms of seasonal affective disorder (SAD) appear.

As a prevention or treatment for SAD, it is advised to carry out sessions in the morning after waking up.

It is also advised to do one session every morning in 3-week cycles during the dark periods of the year. These figures are to adapt to each person depending on the results after a few days of treatment.

White light maximum intensity option: 30 mins at 30 cm / 45 mins at 50-70cm

Blue light maximum intensity option: 45 mins at 30 cm / 60 mins at 50 cm

Generally, all people wanting to compensate their lack of light (SAD, working in a confined room...) with light therapy must give priority to morning exposure: we recommend an exposure of 30 minutes to 2 hours (see above for the distance) at maximal intensity. This morning exposure ideally just after waking up is recommended for people with delayed sleep phase syndrome (falling asleep late and difficulties waking up in the morning) and for people preparing a trip easterly 2 or 3 days before departure (crossing more than 3 time zones), the effects of jet-lag will be notably minimized.

Evening light therapy is reserved for people working at night, people with advanced sleep phase syndrome (falling asleep too early and waking up too early) and for people preparing a trip westerly (crossing more than 3 time zones) to avoid jet-lag.

After your daily light therapy session, you can use your Dayvia Sundesk lamp all day long as a desk

lamp by orienting the illuminating part towards the work surface and reducing the intensity of the white light according to your desires.

## DEVICE MAINTENANCE

Unplug the lamp and the power supply unit and leave it to cool down before cleaning it.

Never use abrasive substances or detergents to clean the device: use a damp cloth making sure that no liquid penetrates into the device.

Being a LED technology, there is no maintenance and the lifespan is over 30000 hours.

## TECHNICAL CHARACTERISTICS

The characteristics are prone to modifications without notice.

- Power supply: 21W
- Tension of the equipment: 220-240 V
- Tension of the power pack: 6V
- Frequency: 50-60Hz
- Insulation: Class II
- Humidity rate during use: between 15% and 90% (without condensation)
- Penetration of foreign bodies and water: IP21
- Atmospheric pressure during use: between 700 hPa and 1060 hPa
- Storage and transport temperature: between -25 ° C and + 70 ° C with relative humidity up to 90% without condensation and water vapor pressure not exceeding 50 hPa

- Overvoltage category II
- Cord Length: 2m65
- Dimension of product: 32.9 x 36 x 9.5 (W x H x D in cm)
- Weight of the device without power supply: 1.3 kg
- Weight of the device with power supply: 1.5 kg
- Class II A medical device
- CE0459

DAYVIA Sundesk is certified conform to all the guideline norms 93/42/CEE and 2007/47/CEE applied to medical devices in Europe :

- Conform to electrical security
- Conform to photobiological risks of light : group 0 (no risk)
- Conform to CEM (magnetic fields) : Class B and Group 1 (Emission RF CISPR 11). Electrostatic discharge immunity test: 15kV for air discharge, 8 Kv for contact discharge. Radiated, radio-frequency, electromagnetic field immunity test: 10V/m. Electrical fast transient/burst immunity test: 2 kV. Surge immunity test:  $\pm 1$  kV differential mode and  $\pm 2$  kV common mode. Immunity to conducted disturbances, induced by radio-frequency fields: 150kHz to 80MHz 3V, 6V for other frequency bands. Power frequency magnetic field immunity test: 3A/m. Voltage dips, short interruptions and voltage variations immunity tests: 0% 10ms, 70% 500ms, 0% 5000ms.

# GUARANTEE AND SUPPORT

If you need information or support, please visit the DAYVIA website at <http://www.dayvia.com/>.

## EXPLANATION OF SYMBOLS

Consult the supplied user's guide.

Manufacturer.

Catalogue number

Batch identification.

This device is class II.

This device must not be discarded with household wastes, take it to a waste disposal facility where it can be recycled. By doing so, you will contribute to the protection of the environment.

Equipment with safety insulation.

**IP21**

Degree of protection against liquid and particles.


Temperature limits.

Humidity limits.

This device is certified conform to the guideline regulating medical devices.

**DAYVIA®**  
LUMINOTHERAPIE

Immeuble Le leeds  
253 Bd du Leeds  
59777 LILLE

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*D – Notice DAYVIA Sundesk - Rev B - novembre 2018*