INNOLUX

BRIGHT

VALOVOIMA

Käyttäjän opas
Bruksanvisningar
User's Instructions
Gebrauchsanleitung
Manuel de l'utilisateur
Manuale utente
Gebruiksaanwijzing
Brugervejledning
Brukerveiledning
Kasutusjuhend
Vartotojo vadovas
Инструкция для потребителей
Lietošanas instrukcija
Pomoc użytkownika

Suomi.....	4
Ruotsi.....	6
English.....	8
Deutsch.....	10
Français.....	12
Italiano.....	14
Nederlands.....	16
Dansk.....	18
Norsk.....	20
Eesti keeles.....	22
Polskie.....	24
Русский.....	26
Latviešu	28
Lietuvių k.	30


INNO LUX BRIGHT

INNOJOK Oy
Sirrikuja 3 L
00940 Helsinki
Finland
+ 358 9 4789 2200
www.innolux.fi



Valovoima


1


2



3

4

5

MEDICAL DEVICE*

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WHAT IS BRIGHT LIGHT THERAPY?

Congratulations for purchasing INNOLUX bright light therapy lamp. Before using the lamp, it is good to know what effects bright light has on your body and how to use bright light therapy effectively. Many people suffer from lethargy and depression during the dark winter season. This condition is known as SAD (Seasonal Affective Disorder). Scientific research has shown that as the amount of natural light decreases, the body may physiologically respond to these “unfavourable living conditions” by adjusting certain hormone levels and biological activities. SUB SAD or the Winter Blues is a milder form of SAD affecting a considerably larger number of people. The effects of these adjustments may include:

- Intense fatigue, exhaustion and various sleep disorders
- Increased appetite leading to weight gain
- General passivity and depression
- Decreased work efficiency and loss of vitality

On a sunny day the strength of light is approximately 50 000 – 100 000 lux, whereas indoors it is as little as 100-200 lux. An effective bright light therapy session requires the eyes to be exposed to at least 2 500 lux. Via the retina and the optic nerves the light affects the brain's light sensitive region (suprachiasmatic cellkernel), which adjusts the body's physiological cycles in accordance with the rhythm of day and night. The light reduces the amount of the sleep hormone melatonin, secreted by the pineal gland. Melatonin is a key factor in the regulation of alertness and sleep cycles, but also in the prevention of SAD. Bright light therapy helps you alleviate these problems and regain vitality. If the symptoms are severe, they should be treated under medical supervision.

THERAPY INSTRUCTIONS

Start bright light therapy before any symptoms of SAD appear, preferably soon after nights start drawing in. Take a bright light session 3-5 times a week throughout the winter. If the symptoms are not relieved or they reappear, try taking daily sessions for 1-2 weeks. You may even begin with daily sessions or repeat this intensity whenever you feel you need it.

One session lasts for 0,5-2 hours. Longer than two hours does not intensify the effect but there's no harm of longer sessions. It is advisable to take the session between 6.00 a.m. and 6.00 p.m. but it is most effective in the morning between 6 and 8. A session in the evening may disturb sleep. If you tend to have evening sleepiness (common to elderly people), you might want to try taking the sessions earlier in the evening in order to stay up longer. You can also divide the session into several shorter periods during the same day if it suits your schedule better.

There is no need to look straight into the light source, but you should keep it within your visual field during the whole session. You can read, write or do other things in the meantime. The effectiveness of bright light therapy correlates with the distance between your face and the light source. Make sure that your face is not further from the lamp than the announced therapy distance on the model instructions. You can shorten the required session time by sitting closer to the lamp.

Bright light therapy also helps speed recovery from jet-lag. A session in the evening will delay your sleep cycle, a session in the morning will have an opposite effect.

IMPORTANT NOTES ON BRIGHT LIGHT THERAPY

The lamp comes with an electronic ballast which raises the light frequency up to 30 000 Hz, where flicker is not noticeable to the eye. After the first few sessions, mild headache or eye strain may occur. This is quite normal and such side effects usually disappear in a few days. If this does not happen, stop the therapy and consult your doctor before continuing. Not all people suffer from SAD due to lack of daylight. Exhaustion, sleep disorders and depression may have other causes. Bright light therapy is most effective when combined with regular physical exercise and a healthy diet.

Consult your doctor before starting bright light therapy if you suffer from

- eye illnesses (e.g. eye inflammation or hypersensitivity)
- mental disorders
- or if you use psychotropic drugs

When you start the bright light therapy, it is useful to take notes on the effects it has on you. This helps you determine which way the therapy works best for you. The need for bright light therapy as well as its duration and frequency vary individually. Use the therapy diary on the back of this booklet.

GETTING STARTED

1. Follow the device-specific instructions on page 3. The LED light sources are in place, and do not need installation.
2. The bright light device is designed for dry premises and must not be placed e.g. in a bathroom, sauna, or other wet or high-humidity area. Similarly, the device should never be used outdoors.
3. Plug the power cord into the socket. The voltage must be 200-240 V.
4. Turn the power on using the switch on the connecting wire.
5. Turn the device off after use.

LED LIGHT SOURCES

In addition to being energy-efficient, LED light sources have a long lifespan. Unlike traditional light sources, LED light sources do not need replacing. LED light source panel is ready installed in the light fixture. The diffuser must not be opened by the user. If the diffuser comes off for any reason, take care not to touch the LED light sources, as this may damage them.

TROUBLESHOOTING

1. If the light fails to work, ensure that the switch is turned on. If the light still won't turn on, contact the retailer or manufacturer for help.
2. LED light sources must be replaced by professional electricians only. If the light sources have been damaged, contact a professional electrician, retailer, or manufacturer

SERVICING THE DEVICE

Keep the device free of dust and dirt. Ensure that the device is not plugged in to the mains, and that it has cooled down after use. When required, wipe all exterior surfaces with a lightly moistened cloth. Do not use abrasive agents. The device must be kept away from water.

IMPORTANT

- Do not let children play with the device.
- Damaged power cord must be replaced by a professional electrician only.
- LED light sources must be replaced by a professional electrician only

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